

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



ANATOMICAL FIXATION SOLUTIONS

Small fragment | Large fragment | Cable plate | Osteotomy
Pediatric | Mini plates | Screws

CONTENTS

●	About SANOMED	03
●	Vision & Mission	04
●	SMALL FRAGMENT	05
●	LARGE FRAGMENT	23
●	Cable Plate	31
●	Osteotomy	35
●	Pediatric	38
●	Mini Plates	41
●	SCREWS	45
●	NAILING SYSTEMS	53
●	New Products	65

PROFESSIONAL ORTHOPEDIC PRODUCT GUIDE

Specifications and availability are subject to confirmation.

ABOUT SANOMED

SCREWS PLATES & NAILING

plies provider focused on advanced orthopedic support hospitals and surgeons with and practical technical coordination.

PRODUCT CATALOG 2025

● TRAUMA

Plates, screws and fixation systems

● SPORTS MEDICINE


Selected surgical solutions

● SERVICE

Local supply and professional support

PRECISION. STRENGTH. HEALING.

Orthopedic solutions supported by a responsive local team.



VISION & MISSION

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

Our vision is to improve patient care by connecting clinicians with reliable technologies, consistent service and long-term value.

● MISSION

To deliver carefully selected implants and instruments with dependable logistics, clear product information and responsive support for healthcare professionals.

OUR COMMITMENT


Quality-minded selection | Accurate coordination | Reliable local support

Continuous development | Ethical professional relationships

01

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



SCIAVICI E
ATE
SCREWS
PLATES
& NAILING


*INDICATION:
Shaft fractures of
the clavicle bone

PRODUCT CATALOG 2025


K- wire holes on both ends to aid in
plate positioning.

Anatomical design
matches the S-shape of
the clavicle

Combi holes for
both locking & non
locking screws




Holes	Length (mm)	Right	Left
4	68	D02T.111R104	D02T.111L104
6	95	D02T.111R106	D02T.111L106
8	119	D02T.111R108	D02T.111L108





DISTAL CLAVICLE SCREWS PLATES & NAILING

*INDICATION:
Lateral third clavicle
fractures

PRODUCT CATALOG 2025



Multiple distal screw configurations to achieve construct stability



Holes	Length (mm)	Right	Left
3	68	D02T:117R103	D02T:117L103
4	82	D02T:117R104	D02T:117L104
5	96	D02T:117R105	D02T:117L105
6	110	D02T:117R106	D02T:117L106
7	124	D02T:117R107	D02T:117L107
8	138	D02T:117R108	D02T:117L108

ACROMIOCLAVICULATE SCREWS PLATES & NAILING

*INDICATION:
Fixation for lateral
clavicle fractures &
acromioclavicular
joint injuries

mm K- wire hole

PRODUCT CATALOG 2025

Combi holes for
both locking & Non
locking screws


Hook provide additional; support to both
lateral clavicle fractures & acromioclavicular
joint dislocations.



Smooth hook design & rounded
plate edges that minimizes soft
tissue irritaions



3.5 mm cortical locking screw



3.5 mm cortical non locking screw

Holes	Length (mm)	Right	Left
3	66	D02T.116R103	D02T.116L103
4	79	D02T.116R104	D02T.116L104
5	91	D02T.116R105	D02T.116L105
6	104	D02T.116R106	D02T.116L106
7	117	D02T.116R107	D02T.116L107
8	129	D02T.116R108	D02T.116L108



PROXIMAL HUMERUS

*INDICATION:
Complex fracture,
fracture
dislocations, osteotomies & nonunions of
the proximal
humerus

SCREWS PLATES & NAILING


PRODUCT CATALOG 2025

Ten suture holes around the
perimeter of the proximal end
including undercuts that facilitates
plate positioning

Elongated combi hole
that helps in plate
positioning



Plate head allows to
insert 9 screws 4.0 mm
cancellous locking



-  3.5 mm cortical locking screw
-  3.5 mm cortical non locking screw
-  4.0 mm cancellous locking screw

Holes	Length (mm)	Universal
3	90	D02T.138103
4	102	D02T.138104
5	114	D02T.138105
6	126	D02T.138106
7	138	D02T.138107
8	150	D02T.138108
9	162	D02T.138109
10	174	D02T.138110
11	186	D02T.138111
12	198	D02T.138112
13	210	D02T.138113

DOIT Humerus Plate SCREWS PLATES & NAILING

*INDICATION:
Temporary Fixation,
correction or
stabilization
of Humerus




PRODUCT CATALOG 2025




Undercuts of the
LC-DCP limit contact
between plate & bone

Combi holes for both
locking & Non locking
screws



A mid-section which is placed
at the position of fragments to
be compressed



3.5 mm cortical locking screw



3.5 mm cortical non locking screw

Holes	Length (mm)	Universal
6	118	DOIT.104106
7	136	DOIT.104107
8	154	DOIT.104108
9	172	DOIT.104109
10	190	DOIT.104110

DISTAL LATERAL

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025


allow locking or compression options

ATE

*INDICATION:
Intra-articular &
supracondylar
fractures of the
distal humerus



Combi holes for both
locking & Non locking
screws



Undercuts limit
plate-to-bone contact to
reduce impairment of
blood supply



Pre-contoured design for anatomic
fit.
No need to bend during the
surgery



2.7 mm cortical locking screw




3.5 mm cortical locking screw



3.5 mm cortical non locking screw

Holes	Length (mm)	Right	Left
3	64	D02T.118R103	D02T.118L103
4	77	D02T.118R104	D02T.118L104
5	90	D02T.118R105	D02T.118L105
6	103	D02T.118R106	D02T.118L106
7	116	D02T.118R107	D02T.118L107
8	129	D02T.118R108	D02T.118L108
9	142	D02T.118R109	D02T.118L109



DISTAL MEDIAL

SCREWS PLATES & NAILING

ATE

*INDICATION:
Intra-articular &
supracondylar
fractures of the
distal humerus

PRODUCT CATALOG 2025



*Pre-contoured design for anatomic fit.
No need to bend during the surgery



3.5 mm cortical locking screw



3.5 mm cortical non locking screw



Elongated combi hole in
the neck and shaft
facilitates plate adjustment
and allows locking or
compression options



*Undercuts limit
plate-to-bone contact to
reduce impairment of
blood supply

Holes	Length (mm)	Right	Left
3	58	D02T.122R103	D02T.122L103
4	71	D02T.122R104	D02T.122L104
5	84	D02T.122R105	D02T.122L105
6	97	D02T.122R106	D02T.122L106
7	110	D02T.122R107	D02T.122L107
8	123	D02T.122R108	D02T.122L108
9	136	D02T.122R109	D02T.122L109




PROXIMAL RADIUS SCREWS PLATES & NAILING

INDICATION: Extra-articular and intra-articular fractures of the proximal radius and multi-fragmented radial neck fractures.

PRODUCT CATALOG 2025

Indication for
distraction



Combi holes allow fixation with locking screws in the threaded section for angular stability, and cortex screws in the Dynamic Compression Unit (DCU) section for distraction



2.4 - 2.7 mm cortical locking screw



2.4 - 2.7 mm Non-locking screw

Holes	Length (mm)	Right	Left
2	36	D02T.151R02	D02T.151L02
3	44	D02T.151R03	D02T.151L03
4	53	D02T.151R04	D02T.151L04


OLECRANON ATE SCREWS PLATES & NAILING

*INDICATION:

Complex extra- and intra-articular olecranon fractures, pseudoarthrosis of the proximal ulna and simple olecranon fractures

PRODUCT CATALOG 2025

facilitates plate adjustment and allow locking or compression options




The proximal spoon-shaped part of the plate is slightly thinner than the shaft




Combi holes for both locking & Non locking screws



*Pre-contured design for anatomic fit. No need to bend during the surgery



Holes	Length (mm)	Right	Left
3	97	D02T.126R103	D02T.126L103
4	110	D02T.126R104	D02T.126L104
5	123	D02T.126R105	D02T.126L105
6	136	D02T.126R106	D02T.126L106
7	149	D02T.126R107	D02T.126L107
8	162	D02T.126R108	D02T.126L108



ONE-TRIGGER™ ULNAR

*INDICATION: Shaft fractures of Ulna, Radius & Fibula

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



Available with combi
holes allowing fixation
with locking screws and
compression with cortex
screws



*Low plate profile and
rounded plate edges
minimize tendon and soft
tissue irritation




3.5 mm cortical locking screw



3.5 mm cortical non locking screw

Holes	Length (mm)	Universal
4	54	D02T.108104
5	66	D02T.108105
6	78	D02T.108106
7	90	D02T.108107
8	102	D02T.108108
9	114	D02T.108109
10	126	D02T.108110
11	138	D02T.108111
12	150	D02T.108112



LC DCP II NA RADIUS TE

SCREWS PLATES & NAILING

*INDICATION:
Temporary Fixation,
correction or
stabilization
of Ulna & Radius

PRODUCT CATALOG 2025


screws



2-mm K-wire hole



*Undercuts of the
LC-DCP limit contact
between plate & bone



A mid-section which is placed
at the position of fragments to
be compressed




3.5 mm cortical locking screw



3.5 mm cortical non locking screw

Holes	Length (mm)	Universal
4	61	D02T.102104
5	74	D02T.102105
6	87	D02T.102106
7	100	D02T.102107
8	113	D02T.102108
9	126	D02T.102109
10	139	D02T.102110
12	165	D02T.102112
14	191	D02T.102114



DISTAL RADIUS



SCREWS
PLATES
& NAILING

PRODUCT CATALOG 2025

3 PLATE

INDICATION: Fixation of complex intra- and extra-articular fractures and osteotomies of the distal radius and other small bones

The window allows
resection of fracture lines
to allow fine contouring.



Holes	Length (mm)	Right	Left
3	57	D02T.150R03	D02T.150L03
4	65	D02T.150R04	D02T.150L04
5	73	D02T.150R05	D02T.150L05

Holes	Length (mm)	Right	Left
3	57	D02T.153R03	D02T.153R03
4	65	D02T.153R04	D02T.153R04
5	73	D02T.153R05	D02T.153R05



DISTAL MEDIAL SCREWS PLATES & NAILING

***INDICATION:**
Fixation of intra- and
extra-articular
fractures
and osteotomies of
the distal tibia

PRODUCT CATALOG 2025

Elongated holes in shaft
aid in initial plate
positioning





Distal tab for optional
medial malleolus screw
accepts 3.5 mm locking,
3.5 mm cortex, 4.0 mm
cortex



Plates can be used with the
LISS Plate technique and
aiming instrumentation

*The plate Available With shaft
combi holes that accepts 5.0
cortical locking & 4.5 cortical
non locking screws

Holes	Length (mm)	Right	Left
4	126	D01T.164R04	D01T.164L04
6	158	D01T.164R06	D01T.164L06
8	190	D01T.164R08	D01T.164L08
10	222	D01T.164R10	D01T.164L10
12	254	D01T.164R12	D01T.164L12
14	286	D01T.164R14	D01T.164L14



DISTAL LATERAL KING

*INDICATION: Extra-articular and intra-articular distal tibia fractures

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025


the head and 3.5mm screws in the shaft



*Three types are available



*Four K-wire holes in the head accept 2.0 mm Kirschner wires



Accepts 3.5mm screws in the head and the shaft



2.7 mm cortical locking screw




3.5 mm cortical locking screw



3.5 mm cortical non locking screw

Holes	Length (mm)	Right	Left
3	74	D01T.132R103	D01T.132L103
4	87	D01T.132R104	D01T.132L104
5	100	D01T.132R105	D01T.132L105
6	113	D01T.132R106	D01T.132L106
7	126	D01T.132R107	D01T.132L107
9	152	D01T.132R109	D01T.132L109




CALCANEUS

SCREWS
PLATES
& NAILING

PRODUCT CATALOG 2025




Large number of holes
address multiple
fracture patterns



PLATE

INDICATION: Fractures and osteotomies of the calcaneus including, extra-articular, intra-articular, joint depression and severely comminuted fractures

omical
ponding to
r of the*It is available in three
es (Small,
ge)



It is available in three
different sizes (Small,
Medium, Large)



3.5 mm cortical locking screw



4.0 mm cancellous locking screw

Holes	Length (mm)	Right	Left
Small	60	D01T.161RS	D01T.161LS
Medium	74	D01T.161RM	D01T.161LM
Large	80	D01T.161RL	D01T.161LL


CURVED

SCREWS PLATES & NAILING

TION TE

*INDICATION:
Temporary fixation,
correction or
stabilization of
bones in the pelvis

PRODUCT CATALOG 2025



3.5 mm cortical locking screw



4.0 mm cancellous locking screw

Accepts 3.5mm cortical
locking screws or 4.0
cancellous locking screws

Holes	Length (mm)	Universal
4	59	D07T.109104
5	71	D07T.109105
6	82	D07T.109106
7	93	D07T.109107
8	104	D07T.109108
9	115	D07T.109109
10	125	D07T.109110
11	135	D07T.109111
12	144	D07T.109112
13	153	D07T.109113
14	161	D07T.109114



STRAIGHT SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

***INDICATION:**
Temporary fixation,
correction or
stabilization of
bones in the pelvis

Low profile and rounded
edges minimize soft tissue
irritation

Combi holes for both
locking & non locking
screws



A mid-section which is
placed at the position of
fragments to be
compressed



3.5 mm cortical locking screw



3.5 mm cortical non locking screw



4.0 mm cancellous locking screw

Holes	Length (mm)	Universal
4	57	D01T.17604
5	72	D01T.17605
6	86	D01T.17606
7	101	D01T.17607
8	115	D01T.17608
9	130	D01T.17609
10	144	D01T.17610
11	159	D01T.17611
12	173	D01T.17612
13	188	D01T.17613
14	202	D01T.17614

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



PROXIMAL FEMUR ATE

*INDICATION: Repairing multifragmentary fractures, malunions, periprosthetic fractures and osteotomies of proximal femur


SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

mm cannulated locking screw





*Undercuts of the LC-DCP limit contact between plate & bone



Anatomically contoured to approximate the lateral aspect of the proximal femur


Elongated Cambi hole in the neck and shaft facilitate plate adjustment and allow locking or compression options



6.5 mm cannulated cancellous locking screw



5.0 mm cortical locking screw



4.5 mm cortical non locking screw

Holes	Length (mm)	Right	Left
3	120	D01T:171R03	D01T:171L03
4	138	D01T:171R04	D01T:171L04
5	156	D01T:171R05	D01T:171L05
6	174	D01T:171R06	D01T:171L06
7	192	D01T:171R07	D01T:171L07
8	210	D01T:171R08	D01T:171L08
9	228	D01T:171R09	D01T:171L09
10	246	D01T:171R10	D01T:171L10
11	264	D01T:171R11	D01T:171L11
12	282	D01T:171R12	D01T:171L12

DISTAL FEMUR ATE

SCREWS PLATES & NAILING

*INDICATION:
Fractures in normal
or osteopenic bone;
nonunions
and malunions; and
osteotomies of the
femur

PRODUCT CATALOG 2025

Kirschner wire holes
accept 2.0 mm Kirschner
wires to temporarily fix the
plate

Combi holes for both
locking & non locking
screws



Distal head accepts 5.5
mm cancellous locking
screw



5.0 mm cortical locking screw



4.5 mm cortical non locking screw



5.5 mm cancellous locking screw

Holes	Length (mm)	Right	Left
5	192	D01T:144R105	D01T:144L105
6	217	D01T:144R106	D01T:144L106
7	242	D01T:144R107	D01T:144L107
8	267	D01T:144R108	D01T:144L108
9	297	D01T:144R109	D01T:144L109
10	317	D01T:144R110	D01T:144L110
11	342	D01T:144R111	D01T:144L111
12	367	D01T:144R112	D01T:144L112
13	392	D01T:144R113	D01T:144L113



LC DCP BROAD FEMUR ATE

INDICATION: Fixation
of various long bones,
such as the humerus
& femur


SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

A mid-section which is
placed at the position of
fragments to be
compressed



Undercuts of the
LC-DCP limit contact
between plate & bone



5.0 mm cortical locking screw



4.5 mm cortical non locking screw

Holes	Length (mm)	Universal
4	85	D01T.102104
5	103	D01T.102105
6	121	D01T.102106
7	139	D01T.102107
8	157	D01T.102108
9	175	D01T.102109
10	193	D01T.102110
11	211	D01T.102111
12	229	D01T.102112
13	247	D01T.102113
14	265	D01T.102114
15	283	D01T.102115



PROXIMAL LATERAL TIBIA PLATE

*INDICATION:
Stabilization of fractures
of the proximal tibia

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



Plates are pre-contoured
to create a fit that requires
little or no additional
bending



Three 2.0 mm holes for
preliminary fixation
with Kirschner wires, or
meniscal repair with
sutures



5.0 mm cortical locking screw



4.5 mm cortical non locking screw



5.5 mm cancellous locking screw

Holes	Length (mm)	Right	Left
5	141	D01T:157R105	D01T:157L105
7	175	D01T:157R107	D01T:157L107
9	209	D01T:157R109	D01T:157L109
11	243	D01T:157R111	D01T:157L111
13	277	D01T:157R113	D01T:157L113
15	311	D01T:157R115	D01T:157L115

LC-DCP NARROW


SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

*INDICATION: Fixation of various long bones, such as the tibia



Undercuts of the LC-DCP limit contact between plate & bone



Amid-section which is placed at the position of fragments to be compressed



5.0 mm cortical locking screw



4.5 mm cortical non locking screw

Holes	Length (mm)	Universal
4	85	D01T:107104
5	103	D01T:107105
6	121	D01T:107106
7	139	D01T:107107
8	157	D01T:107108
9	175	D01T:107109
10	193	D01T:107110
11	211	D01T:107111
12	229	D01T:107112
13	247	D01T:107113
14	265	D01T:107114
15	283	D01T:107115
16	301	D01T:107116

TRI ITTDECC ATE

SCREWS PLATES & NAILING

*INDICATION: Used to rigidly hold in place fractures at the end of long bones, especially at the knee

PRODUCT CATALOG 2025

bending



Elongated Cambi hole in the neck and shaft facilitate plate adjustment and allow locking or compression options

*It is Ideal to buttress small metaphyseal fragment with the diaphysis



5.0 mm cortical locking screw



4.5 mm cortical non locking screw



5.5 mm cancellous locking screw

Holes	Length (mm)	Universal
3	90	D01T:140103
4	106	D01T:140104
5	122	D01T:140105
6	138	D01T:140106
7	154	D01T:140107
8	170	D01T:140108
9	186	D01T:140109

IRITTDECC

SCREWS
PLATES
& NAILING

PRODUCT CATALOG 2025

Longer diameter in the neck and shaft facilitate plate adjustment and allow locking or compression options



*It is Ideal to buttress small metaphyseal fragment with the diaphysis



Plate is broadened and carefully contoured at the joint end of the plate.



5.0 mm cortical locking screw



4.5 mm cortical non locking screw



5.5 mm cancellous locking screw


Holes	Length (mm)	Right	Left
3	76	D01T.141R103	D01T.141L103
4	92	D01T.141R104	D01T.141L104
5	108	D01T.141R105	D01T.141L105
6	124	D01T.141R106	D01T.141L106
7	140	D01T.141R107	D01T.141L107
8	156	D01T.141R108	D01T.141L108

03

TE

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



PROXIMAL FEMUR

SCREWS
PLATES
& NAILING

PRODUCT CATALOG 2025

Anatomically pre-contoured to approximate the lateral aspect of the proximal femur



45mm Trochanteric Grip Plate

ATE

*INDICATION: The complex fractures, fracture dislocations, osteotomies and nonunions of the femur



Two proximal hooks engage the superior tip of the greater trochanter



Holes on the both sides of the plate accept 2.0 m K-wire



5.0 mm cortical locking screw



4.5 mm cortical non locking screw



5.5 mm cancellous locking screw


Holes	Length (mm)	Right	Left
6	180	D01T.181R180	D01T.181L180
7	213	D01T.181R230	D01T.181L213
8	233	D01T.181R233	D01T.181L233
9	265	D01T.181R265	D01T.181L265
11	285	D01T.181R285	D01T.181L285

DISTAL LATERAL E ATE SCREWS PLATES & NAILING

*INDICATION:
Buttressing multi-
fragmentary distal
femur
fractures

PRODUCT CATALOG 2025


The highest level of compression is provided by this plate through the plate itself and the cable system



Holes on the both sides of the plate accept 2.0 m K-wire.



8 Screw Holes & 3 K-wire holes on the distal head of the plate



Holes	Length (mm)	Right	Left
3	119	D02T.191R03	D02T.191L03
4	139	D02T.191R04	D02T.191L04
5	159	D02T.191R05	D02T.191L05
6	179	D02T.191R06	D02T.191L06
7	199	D02T.191R07	D02T.191L07
8	219	D02T.191R08	D02T.191L08
9	239	D02T.191R09	D02T.191L09

CEMI ID CABLE ATE SCREWS PLATES & NAILING

PRODUCT CATALOG 2025


*INDICATION: The complex fractures, fracture dislocations, osteotomies and nonunions of the femur

Holes on the both sides of the plate accept 2.0 m K-wire



*The highest level of compression is provided by this plate through the plate itself and the cable system

Holes	^(mm) Length	Universal
5	133	D02T.191R05
6	158	D01T.178106
7	183	D01T.178107
8	208	D01T.178108
9	233	D01T.178109
11	283	D01T.178111



04

Y

SCREWS PLATES & NAILING


PRODUCT CATALOG 2025



TIBIAL OSTEOTOMY
SCREWS
PLATES
& NAILING

*INDICATION: Proximal tibia diseases including osteotomies, bone and joint deformities, fractures and malalignment

PRODUCT CATALOG 2025



Without Wedge



With Wedge



Holes on the both sides of the plate accept 2.0 mm K-wire



Pre-contoured design for anatomic fit, minimizing soft tissue irritation



5.0 mm cortical locking screw



4.5 mm cortical locking screw



5.5 mm cancellous locking screw



osteotomy wedge

Holes	Length (mm)	Universal
small	44	D01.1833
medium	48	D01.1834
large	52	D01.1835


CEMI ID OSTEOTOMY

*INDICATION:
Buttressing
multifragmentary distal
femur
fractures

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

Cambi holes allow fixation with locking screws in the threaded section for angular stability, and cortex screws in the Dynamic Compression Unit (DCU) section for distraction



*Pre-contoured design
for anatomic fit,
minimizing soft tissue
irritation



5.0 mm cortical locking screw



4.5 mm cortical locking screw




5.5 mm cancellous locking screw

Holes	Length (mm)	Right	Left
5	142.5	D01T.150R105	D01T.150L105
7	174.5	D01T.150R107	D01T.150R107

05

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



EDIDLVCEAI EIGHT

SCREWS PLATES & NAILING

*INDICATION: Gradual correction of the angled deformities to upper and lower extremities during the adolescence and infant periods

PRODUCT CATALOG 2025




K-wire holes on the center and both end of the plate



Low profile plates

Works with 4.0 & 4.5 mm cannulated cancellous screws



cannulated cancellous non locking screws

4.5 mm cortical non locking screw



Holes	Length (mm)	Universal
12	28	D02T.14012
14	30	D02T.14014
16	32	D02T.14016
18	34	D02T.14018

PEDIATRIC HIP SCREWS PLATES & NAILING


*INDICATION:
Temporary fixation,
correction or
stabilization of
bones in various
anatomical regions in
the infants and children

PRODUCT CATALOG 2025

3.5 mm Pediatric Plates




The proximal parts of the
each type plates designed in
100°-110°-120°-140° for varus
osteotomies, have offset
structure.




5.0 mm Pediatric Plates


Two types of the plates are
available




3.5 mm cortical non locking screw



4.0 mm cancellous locking screw



3.5 mm cortical locking screw



4.5 mm cortical locking screw



5.0 mm cortical locking screw



5.5 mm cancellous locking screw


Holes	Length (mm)	100	110	120	140
3	92	D01T:175L03	D01T:175R03	D01T:175R03	D01T:175M03
4	108	D01T:175L04	D01T:175R04	D01T:175R04	D01T:175M04
5	124	D01T:175L05	D01T:175R05	D01T:175R05	D01T:175M05
6	140	D01T:175L06	D01T:175R06	D01T:175R06	D01T:175M06
7	156	D01T:175L07	D01T:175R07	D01T:175R07	D01T:175M07

Holes	Length (mm)	100	110	120	140
3	73.5	D01T:174L03	D01T:174R03	D01T:174S03	D01T:174M03
4	86.5	D01T:174L04	D01T:174R04	D01T:174S04	D01T:174M04
5	99.5	D01T:174L05	D01T:174R05	D01T:174S05	D01T:174M05
6	112.2	D01T:174L06	D01T:174R06	D01T:174S06	D01T:174M06
7	125.2	D01T:174L07	D01T:174R07	D01T:174S07	D01T:174M07

06

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025




SCREWS
PLATES
& NAILING

PRODUCT CATALOG 2025

ES

*INDICATION: These plates are indicated to address a wide range of small bone fractures including the hand, wrist, foot and ankle




Mesh	H	L	Orbital	Straight	T	Y	X	Z
D03T.115100	D03T.109104	D03T.105103	D03T.113100	D03T.101103	D03T.103103	D03T.107103	D03T.111101	D03T.110101
D03T.115101	D03T.109106	D03T.105104	D03T.113101	D03T.101104	D03T.103104	D03T.107104		
D03T.115102	D03T.109108	D03T.105105	D03T.113102	D03T.101105	D03T.103105	D03T.107105		
		D03T.105106		D03T.101106	D03T.103106	D03T.107106		
		D03T.105107		D03T.101107	D03T.103107	D03T.107107		
		D03T.105108		D03T.101108	D03T.103108	D03T.107108		

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

TES

*INDICATION: These plates are indicated to address a wide range of small bone fractures including the hand, wrist, foot and ankle





Mesh	H	L	Orbital	Straight	T	Y	X	Z
D03T.131102	D03T.125104	D03T.121103	D03T.129100	D03T.117103	D03T.119103	D03T.123103	D03T.129101	D03T.126101
D03T.131101	D03T.125106	D03T.121104	D03T.129101	D03T.117104	D03T.119104	D03T.123104		
D03T.131100	D03T.125108	D03T.121105	D03T.129102	D03T.117105	D03T.119105	D03T.123105		
		D03T.121106		D03T.117106	D03T.119106	D03T.123106		
		D03T.121107		D03T.117107	D03T.119107	D03T.123107		
		D03T.121108		D03T.117108	D03T.119108	D03T.123108		

PLATES

SCREWS
PLATES
& NAILING

PRODUCT CATALOG 2025

*INDICATION: These plates are indicated to address a wide range of small bone fractures including the hand, wrist, foot and ankle




Mesh	H	L	Orbital	Straight	T	Y	X	Z
D03T.131102	D03T.141104	D03T.137103	D03T.145100	D03T.133103	D03T.135103	D03T.139103	D03T.143101	D03T.142101
D03T.131101	D03T.141106	D03T.137104	D03T.145101	D03T.133104	D03T.135104	D03T.139104		
D03T.131100	D03T.141108	D03T.137105	D03T.145102	D03T.133105	D03T.135105	D03T.139105		
		D03T.137106		D03T.133106	D03T.135106	D03T.139106		
		D03T.137107		D03T.133107	D03T.135107	D03T.139107		
		D03T.137108		D03T.133108	D03T.135108	D03T.139108		

07

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025



CORTICAL
LOCKING SCREWS

*Works with star
screwdriver (Torx)

PLATES
& NAILING

PRODUCT CATALOG 2025



1.5 mm
2.0 mm
2.4 mm
2.7 mm
3.5 mm
5.0 mm

1.5 (T5)	2.0 (T6)	2.4 (T8)	2.7 (T8)	3.5 (T15)	5.0 (T25)
D05T.1506	DOST.2006	DOST.2408	DOST.2708	DOST.3S10	DOST.S016
D05T.1507	DOST.2007	DOST.2409	DOST.2710	DOST.3S12	DOST.S018
D05T.1508	DOST.2008	DOST.2410	DOST.2712	DOST.3S14	DOST.S020
D05T.1509	DOST.2009	DOST.2411	DOST.2714	DOST.3S16	DOST.S022
D05T.1510	DOST.2010	DOST.2412	DOST.2716	DOST.3S18	DOST.S024
D05T.1511	DOST.2011	DOST.2414	DOST.2718	DOST.3S20	DOST.S026
D05T.1512	DOST.2012	DOST.2416	DOST.2720	DOST.3S22	DOST.S028
D05T.1513	DOST.2013	DOST.2418	DOST.2722	DOST.3S24	DOST.S030
D05T.1514	DOST.2014	DOST.2420	DOST.2724	DOST.3S26	DOST.S032
D05T.1515	DOST.2015	DOST.2422	DOST.2726	DOST.3S28	DOST.S034
D05T.1516	DOST.2016	DOST.2424	DOST.2728	DOST.3S30	DOST.S036
D05T.1517	DOST.2017	DOST.2426	DOST.2730	DOST.3S32	DOST.S038
D05T.1518	DOST.2018	DOST.2428	DOST.2732	DOST.3S34	DOST.S040
D05T.1519	DOST.2019	DOST.2430	DOST.2734	DOST.3S36	DOST.S042
D05T.1520	DOST.2020	DOST.2432	DOST.2736	DOST.3S38	DOST.S044
	DOST.2022	DOST.2434	DOST.2738	DOST.3S40	DOST.S046
	DOST.2024	DOST.2436	DOST.2740	DOST.3S42	DOST.S048
	DOST.2026	DOST.2438	DOST.2742	DOST.3S44	DOST.S050
	DOST.2028	DOST.2440	DOST.2744	DOST.3S46	DOST.S052
	DOST.2030		DOST.2746	DOST.3S48	DOST.S054
			DOST.2748	DOST.3S50	DOST.S056
			DOST.2750	DOST.3S52	DOST.S058
				DOST.3S54	DOST.S060
				DOST.3S56	DOST.S065
				DOST.3S70	DOST.S070
					DOST.S075
					DOST.S080


CORTICAL NON
LOCKING SCREWS

*Works with hex
screwdriver in 3.5
and 4.5 screws, and
with star screwdriver
(Torx) in 1.5, 2.0,
2.4 and 2.7 screws.

screws.

PLATES
& NAILING

PRODUCT CATALOG 2025



1.5 mm
2.0 mm
2.4 mm
2.7 mm
3.5 mm
5.0 mm

1.5 (T5)	2.0 (T6)	2.4 (T8)	2.7 (T8)	3.5 (A2.5)	4.5 (A3.5)
D05T.CS1506	D05T.CS2006	D05T.CS2408	D05T.CS2706	D05T.CS3510	D05T.CS4514
D05T.CS1507	D05T.CS2007	D05T.CS2409	D05T.CS2707	D05T.CS3512	D05T.CS4516
D05T.CS1508	D05T.CS2008	D05T.CS2410	D05T.CS2708	D05T.CS3514	D05T.CS4518
D05T.CS1509	D05T.CS2009	D05T.CS2411	D05T.CS2709	D05T.CS3516	D05T.CS4520
D05T.CS1510	D05T.CS2010	D05T.CS2412	D05T.CS2710	D05T.CS3518	D05T.CS4522
D05T.CS1511	D05T.CS2011	D05T.CS2414	D05T.CS2712	D05T.CS3520	D05T.CS4524
D05T.CS1512	D05T.CS2012	D05T.CS2416	D05T.CS2714	D05T.CS3522	D05T.CS4526
D05T.CS1513	D05T.CS2013	D05T.CS2418	D05T.CS2716	D05T.CS3524	D05T.CS4528
D05T.CS1514	D05T.CS2014	D05T.CS2420	D05T.CS2718	D05T.CS3526	D05T.CS4530
D05T.CS1516	D05T.CS2015	D05T.CS2422	D05T.CS2720	D05T.CS3528	D05T.CS4532
D05T.CS1518	D05T.CS2016	D05T.CS2424	D05T.CS2722	D05T.CS3530	D05T.CS4534
D05T.CS1520	D05T.CS2017	D05T.CS2426	D05T.CS2724	D05T.CS3532	D05T.CS4536
D05T.CS1522	D05T.CS2018	D05T.CS2428	D05T.CS2726	D05T.CS3534	D05T.CS4538
D05T.CS1524	D05T.CS2019	D05T.CS2430	D05T.CS2728	D05T.CS3536	D05T.CS4540
D05T.CS1526	D05T.CS2020	D05T.CS2432	D05T.CS2730	D05T.CS3538	D05T.CS4542
D05T.CS1528	D05T.CS2022	D05T.CS2434	D05T.CS2732	D05T.CS3540	D05T.CS4544
D05T.CS1530	D05T.CS2024	D05T.CS2436	D05T.CS2734	D05T.CS3542	D05T.CS4546
	D05T.CS2026	D05T.CS2438	D05T.CS2736	D05T.CS3544	D05T.CS4548
	D05T.CS2028	D05T.CS2440	D05T.CS2738	D05T.CS3546	D05T.CS4550
	D05T.CS2030		D05T.CS2740	D05T.CS3548	D05T.CS4552
			D05T.CS2748	D05T.CS3550	D05T.CS4554
			D05T.CS2750	D05T.CS3552	D05T.CS4556
				D05T.CS3554	D05T.CS4558
				D05T.CS3556	D05T.CS4560
				D05T.CS3558	D05T.CS4565
				D05T.CS3560	D05T.CS4570
				D05T.CS3565	D05T.CS4575
					D05T.CS4580

CANCELLOUS
LOCKING SCREWS

*Works with star
screwdriver (Torx)

PLATES
& NAILING

PRODUCT CATALOG 2025



4.0 mm
5.5 mm


4.0 (T15)	5.5 (T25)
D05T.SL4010	D05T.SL5020
D05T.SL4012	D05T.SL5022
D05T.SL4014	D05T.SL5024
D05T.SL4016	D05T.SL5026
D05T.SL4018	D05T.SL5028
D05T.SL4020	D05T.SL5030
D05T.SL4022	D05T.SL5032
D05T.SL4024	D05T.SL5034
D05T.SL4026	D05T.SL5036
D05T.SL4028	D05T.SL5038
D05T.SL4030	D05T.SL5040
D05T.SL4035	D05T.SL5042
D05T.SL4040	D05T.SL5044
D05T.SL4045	D05T.SL5046
D05T.SL4050	D05T.SL5048
D05T.SL4055	D05T.SL5050
D05T.SL4060	D05T.SL5055
D05T.SL4065	D05T.SL5060
D05T.SL4070	D05T.SL5065
D05T.SL4075	D05T.SL5070
	D05T.SL5075
	D05T.SL5080
	D05T.SL5085
	D05T.SL5090

CANCELLOUS NON
LOCKING SCREWS

*Works with hex
screwdriver

PLATES
& NAILING

PRODUCT CATALOG 2025



4.5 Half	4.5 Full	6.5 Half	6.5 Full
D05T.SH4030	D05T.S4010	D05T.SH65030	D05T.S65030
D05T.SH4035	D05T.S4012	D05T.SH65035	D05T.S65035
D05T.SH4040	D05T.S4014	D05T.SH65040	D05T.S65040
D05T.SH4045	D05T.S4016	D05T.SH65045	D05T.S65045
D05T.SH4050	D05T.S4018	D05T.SH65050	D05T.S65050
D05T.SH4055	D05T.S4020	D05T.SH65055	D05T.S65055
D05T.SH4060	D05T.S4022	D05T.SH65060	D05T.S65060
D05T.SH4065	D05T.S4024	D05T.SH65065	D05T.S65065
D05T.SH4070	D05T.S4026	D05T.SH65070	D05T.S65070
	D05T.S4028	D05T.SH65075	D05T.S65075
	D05T.S4030	D05T.SH65080	D05T.S65080
	D05T.S4040	D05T.SH65085	D05T.S65085
	D05T.S4045	D05T.SH65090	D05T.S65090
	D05T.S4050	D05T.SH65095	D05T.S65095
	D05T.S4055	D05T.SH65100	D05T.S65100
	D05T.S4060	D05T.SH65105	D05T.S65105
	D05T.S4065	D05T.SH65110	D05T.S65110
	D05T.S4070	D05T.SH65115	D05T.S65115
		D05T.SH65120	D05T.S65120



CANNULATED CORTICAL
LOCKING SCREWS

*Works with
cannulated star
screwdriver (Torx)

PLATES
& NAILING



PRODUCT CATALOG 2025

2.7 mm
3.5 mm
5.0 mm


2.7 (T8)	3.5 (15)	5.0 (T25)
D05T.CS2710	D05T.C3510	D05T.C5010
D05T.CS2712	D05T.C3512	D05T.C5012
D05T.CS2714	D05T.C3514	D05T.C5014
D05T.CS2716	D05T.C3516	D05T.C5016
D05T.CS2718	D05T.C3518	D05T.C5018
D05T.CS2720	D05T.C3520	D05T.C5020
D05T.CS2722	D05T.C3522	D05T.C5022
D05T.CS2724	D05T.C3524	D05T.C5024
D05T.CS2726	D05T.C3526	D05T.C5026
D05T.CS2728	D05T.C3528	D05T.C5028
D05T.CS2730	D05T.C3530	D05T.C5030
D05T.CS2732	D05T.C3532	D05T.C5032
D05T.CS2734	D05T.C3534	D05T.C5034
D05T.CS2736	D05T.C3536	D05T.C5036
D05T.CS2738	D05T.C3538	D05T.C5038
D05T.CS2740	D05T.C3540	D05T.C5040
D05T.CS2742	D05T.C3542	D05T.C5042
D05T.CS2744	D05T.C3544	D05T.C5044
D05T.CS2746	D05T.C3546	D05T.C5046
D05T.CS2748	D05T.C3548	D05T.C5048
D05T.CS2750	D05T.C3550	D05T.C5050
	D05T.C3552	D05T.C5052
	D05T.C3554	D05T.C5054
	D05T.C3556	D05T.C5056
	D05T.C3558	D05T.C5058
	D05T.C3560	D05T.C5060
	D05T.C3562	D05T.C5062
	D05T.C3564	D05T.C5064
	D05T.C3566	D05T.C5066
	D05T.C3568	D05T.C5068
	D05T.C3570	D05T.C5070

CANNULATED CANCELLOUS
NON LOCKING SCREWS

*Works with
cannulated hex
screwdriver

PLATES
& NAILING

PRODUCT CATALOG 2025



4.5	6.5 Half	6.5 Full	7.3 Half	7.3 Full
D06T.CS4020	D06T.CSH65030	D06T.CS65030	D06T.CSH73050	D06T.CS73045
D06T.CS4022	D06T.CSH65035	D06T.CS65035	D06T.CSH73055	D06T.CS73050
D06T.CS4024	D06T.CSH65040	D06T.CS65040	D06T.CSH73060	D06T.CS73055
D06T.CS4026	D06T.CSH65045	D06T.CS65045	D06T.CSH73065	D06T.CS73060
D06T.CS4028	D06T.CSH65050	D06T.CS65050	D06T.CSH73070	D06T.CS73065
D06T.CS4030	D06T.CSH65055	D06T.CS65055	D06T.CSH73075	D06T.CS73070
D06T.CS4035	D06T.CSH65060	D06T.CS65060	D06T.CSH73080	D06T.CS73075
D06T.CS4040	D06T.CSH65065	D06T.CS65065	D06T.CSH73085	D06T.CS73080
D06T.CS4045	D06T.CSH65070	D06T.CS65070	D06T.CSH73090	D06T.CS73085
D06T.CS4050	D06T.CSH65075	D06T.CS65075	D06T.CSH73095	D06T.CS73090
D06T.CS4055	D06T.CSH65080	D06T.CS65080	D06T.CSH73100	D06T.CS73095
D06T.CS4060	D06T.CSH65085	D06T.CS65085	D06T.CSH73105	D06T.CS73100
	D06T.CSH65090	D06T.CS65090	D06T.CSH73110	D06T.CS73105
	D06T.CSH65095	D06T.CS65095		D06T.CS73110
	D06T.CSH65100	D06T.CS65100		
	D06T.CSH65105	D06T.CS65105		
	D06T.CSH65110	D06T.CS65110		
	D06T.CSH65115	D06T.CS65115		
	D06T.CSH65120	D06T.CS65120		

HEADLESS
COMPRESSION SCREWS

*Works with
cannulated star
screwdriver (Torx)

PLATES
& NAILING



PRODUCT CATALOG 2025

	2.2	2.7	3.5	4.5	6.5
2.2 mm	D04T.HCS2205	D04T.HCS2705	D04T.HCS3510	D04T.HCS4510	D04T.HCS65020
2.7 mm	D04T.HCS2206	D04T.HCS2706	D04T.HCS3512	D04T.HCS4512	D04T.HCS65025
3.5 mm	D04T.HCS2207	D04T.HCS2707	D04T.HCS3514	D04T.HCS4514	D04T.HCS65030
4.5 mm	D04T.HCS2208	D04T.HCS2708	D04T.HCS3516	D04T.HCS4516	D04T.HCS65035
6.5mm	D04T.HCS2209	D04T.HCS2709	D04T.HCS3518	D04T.HCS4518	D04T.HCS65040
	D04T.HCS2210	D04T.HCS2710	D04T.HCS3520	D04T.HCS4520	D04T.HCS65045
	D04T.HCS2211	D04T.HCS2711	D04T.HCS3521	D04T.HCS4522	D04T.HCS65050
	D04T.HCS2212	D04T.HCS2712	D04T.HCS3522	D04T.HCS4524	D04T.HCS65055
	D04T.HCS2213	D04T.HCS2713	D04T.HCS3523	D04T.HCS4526	D04T.HCS65060
	D04T.HCS2214	D04T.HCS2714	D04T.HCS3524	D04T.HCS4528	D04T.HCS65065
	D04T.HCS2215	D04T.HCS2715	D04T.HCS3525	D04T.HCS4530	D04T.HCS65070
	D04T.HCS2216	D04T.HCS2716	D04T.HCS3526	D04T.HCS4531	D04T.HCS65075
	D04T.HCS2217	D04T.HCS2717	D04T.HCS3527	D04T.HCS4533	D04T.HCS65080
	D04T.HCS2218	D04T.HCS2718	D04T.HCS3528	D04T.HCS4535	D04T.HCS65085
	D04T.HCS2219	D04T.HCS2719	D04T.HCS3529	D04T.HCS4537	D04T.HCS65090
	D04T.HCS2220	D04T.HCS2720	D04T.HCS3530	D04T.HCS4539	D04T.HCS65095
	D04T.HCS2221	D04T.HCS2721	D04T.HCS3531	D04T.HCS4540	D04T.HCS65100
		D04T.HCS2722	D04T.HCS3533	D04T.HCS4541	D04T.HCS65105
		D04T.HCS2723	D04T.HCS3535	D04T.HCS4543	D04T.HCS65110
		D04T.HCS2725	D04T.HCS3538	D04T.HCS4545	D04T.HCS65115
		D04T.HCS2727	D04T.HCS3540	D04T.HCS4548	D04T.HCS65120
		D04T.HCS2729	D04T.HCS3543		
		D04T.HCS2731	D04T.HCS3545		
		D04T.HCS2733	D04T.HCS3548		
		D04T.HCS2735	D04T.HCS3550		
		D04T.HCS2740	D04T.HCS3555		
		D04T.HCS2745	D04T.HCS3560		
		D04T.HCS2750			



TRAUMA

NAILING SYSTEMS

Intramedullary fixation portfolio
for upper and lower extremities

Precision. Strength. Healing.





FOCUSED NAILING PORTFOLIO

Stable fixation, clear selection, streamlined surgical planning.

A concise catalogue of SANOMED intramedullary fixation systems covering upper-limb, proximal femur, femoral shaft, tibia, arthrodesis and flexible nailing applications.

01

HUMERUS

Cannulated humeral nailing system

02

PFNA

Proximal femoral nailing system

03

LATERAL FEMUR

Cannulated lateral femoral nailing

04

TIBIA

Cannulated tibial nailing system

05

ARTHRODESIS

Tibial arthrodesis nailing system

06

FLEXIBLE NAILS

Flexible intramedullary nail range

For healthcare professional use. Product availability and regulatory status may vary by market.

HUMERUS CANNULATED NAIL

Code	Size	Length
SN5024207160	7 mm	160
SN5024207185	7 mm	180
SN5024207200	7 mm	200
SN5024207215	7 mm	220
SN5024207245	7 mm	240
SN5024207260	7 mm	260
SN5024207275	7 mm	280
SN5024207305	7 mm	300
SN5024207320	7 mm	320
SN5024208160	8 mm	160
SN5024208185	8 mm	180
SN5024208200	8 mm	200
SN5024208215	8 mm	220
SN5024208245	8 mm	240
SN5024208260	8 mm	260
SN5024208275	8 mm	280
SN5024208305	8 mm	300
SN5024208320	8 mm	320
SN5024209160	9 mm	160
SN5024209185	9 mm	180
SN5024209200	9 mm	200
SN5024209215	9 mm	220
SN5024209245	9 mm	240
SN5024209260	9 mm	260
SN5024209275	9 mm	280
SN5024209305	9 mm	300
SN5024209320	9 mm	320



Humerus Locking Screw

Code	Size	Length
SN5026230020	5 mm	20
SN5026230025	5 mm	25
SN5026230030	5 mm	30
SN5026230035	5 mm	35
SN5026230040	5 mm	40
SN5026230045	5 mm	45
SN5026230050	5 mm	50
SN5026230055	5 mm	55
SN5026230060	5 mm	60
SN5026230065	5 mm	65
SN5026230070	5 mm	70
SN5026230075	5 mm	75
SN5026230080	5 mm	80

Humerus Locking Screw

Code	Size	Length
SN5026235020	3,5 mm	20
SN5026235025	3,5 mm	25
SN5026235030	3,5 mm	30
SN5026235035	3,5 mm	35
SN5026235040	3,5 mm	40
SN5026235045	3,5 mm	45
SN5026235050	3,5 mm	50
SN5026235055	3,5 mm	55
SN5026235060	3,5 mm	60
SN5026235065	3,5 mm	65
SN5026235070	3,5 mm	70
SN5026235075	3,5 mm	75
SN5026235080	3,5 mm	80

Indications

- Stable or unstable fresh fractures, Pathologic fractures, Refractures, delayed unions and pseudarthrosis.

Specifications

- Humerus Nail Instrument Set

Humerus End Cap

Code	Size
SN5027030000	0
SN5027030005	5 mm
SN5027030010	10 mm



PFNA NAIL

Code	Size	Side	Length
------	------	------	--------

SN5064109170	9 mm	-	180
SN5064109200	9 mm	-	200
SN5064109240	9 mm	-	240
SN5064110170	10 mm	-	180
SN5064110200	10 mm	-	200
SN5064110240	10 mm	-	240
SN5064111170	11 mm	-	180
SN5064111200	11 mm	-	200
SN5064111240	11 mm	-	240
SN5064120180	12 mm	-	180
SN5064120200	12 mm	-	200
SN5064120240	12 mm	-	240
SN5064113170	13 mm	-	180
SN5064113200	13 mm	-	200
SN5064113240	13 mm	-	240

SN5064110300	10 mm	R	300
SN5064110320	10 mm	R	320
SN5064110340	10 mm	R	340
SN5064110360	10 mm	R	360
SN5064110380	10 mm	R	380
SN5064110400	10 mm	R	400
SN5064210300	10 mm	L	300
SN5064210320	10 mm	L	320
SN5064210340	10 mm	L	340
SN5064210360	10 mm	L	360
SN5064210380	10 mm	L	380
SN5064210400	10 mm	L	400

SN5064211300	11 mm	R	300
SN5064211320	11 mm	R	320
SN5064211340	11 mm	R	340
SN5064211360	11 mm	R	360
SN5064211380	11 mm	R	380
SN5064211400	11 mm	R	400
SN5064311300	11 mm	L	300
SN5064311320	11 mm	L	320
SN5064311340	11 mm	L	340
SN5064311360	11 mm	L	360
SN5064311380	11 mm	L	380
SN5064311400	11 mm	L	400



Code	Size	Side	Length
------	------	------	--------

SN5064112300	12 mm	R	300
SN5064112320	12 mm	R	320
SN5064112340	12 mm	R	340
SN5064112360	12 mm	R	360
SN5064112380	12 mm	R	380
SN5064112400	12 mm	R	400
SN5064212300	12 mm	L	300
SN5064212320	12 mm	L	320
SN5064212340	12 mm	L	340
SN5064212360	12 mm	L	360
SN5064212380	12 mm	L	380
SN5064212400	12 mm	L	400

SN5064313300	13 mm	R	300
SN5064313320	13 mm	R	320
SN5064313340	13 mm	R	340
SN5064313360	13 mm	R	360
SN5064313380	13 mm	R	380
SN5064313400	13 mm	R	400
SN5064213300	13 mm	L	300
SN5064213320	13 mm	L	320
SN5064213340	13 mm	L	340
SN5064213360	13 mm	L	360
SN5064213380	13 mm	L	380
SN5064213400	13 mm	L	400

Specifications

- Short and Long PFNA options
- Proximal Femur Nail Instrument Set
- Additional Set for Long PFNA
- Available 2 Compression Lag Screw Options

Indications

- Pertrochanteric fractures, Intertrochanteric fractures, High subtrochanteric fractures, Low subtrochanteric fractures, Femoral shaft fractures, Isolated or combined medial femoral neck fractures.

PFNA NAIL

PFNA Compression Locking Screw

Code	Length
SN5066030065	65
SN5066030070	70
SN5066030075	75
SN5066030080	80
SN5066030085	85
SN5066030090	90
SN5066030095	95
SN5066030100	100
SN5066030105	105
SN5066030110	110
SN5066030115	115
SN5066030120	120



PFNA Cannulated Lag Screw

Code	Length
SN5066140065	65
SN5066140070	70
SN5066140075	75
SN5066140080	80
SN5066140085	85
SN5066140090	90
SN5066140095	95
SN5066140100	100
SN5066140105	105
SN5066140110	110
SN5066140115	115
SN5066140120	120



PFNA Anti-Rotation Locking Screw

Code	Length
SN7938140065	65
SN7938140070	70
SN7938140075	75
SN7938140080	80
SN7938140085	85
SN7938140090	90
SN7938140095	95
SN7938140100	100
SN7938140105	105
SN7938140110	110
SN7938140115	115
SN7938140120	120



PFNA NAIL

PFNA Locking Screw

Code	Size	Length
SN5006230020	5 mm	20
SN5006230025	5 mm	25
SN5006230030	5 mm	30
SN5006230035	5 mm	35
SN5006230040	5 mm	40
SN5006230045	5 mm	45
SN5006230050	5 mm	50
SN5006230055	5 mm	55
SN5006230060	5 mm	60
SN5006230065	5 mm	65
SN5006230070	5 mm	70
SN5006230075	5 mm	75
SN5006230080	5 mm	80
SN5006230085	5 mm	85
SN5006230090	5 mm	90
SN5006230095	5 mm	95
SN5006230100	5 mm	100
SN5006230105	5 mm	105
SN5006230110	5 mm	110
SN5006230115	5 mm	115
SN5006230120	5 mm	120



PFNA Nail End Cap

Code	Size
SN5007030000	0
SN5007030005	5 mm
SN5007030010	10 mm
SN5007030015	15 mm
SN5007030020	20 mm



PFNA Nail End Cap Insert Locking Screw

Code	Size
SN5007030025	-



PFNA Nail End Cap Compression Insert Locking Screw

Code	Size
SN5007030026	-



INTRAMEDULLARY CANNULATED LATERAL FEMUR NAIL

Code	Size	Side	Length		Code	Size	Side	Length
SN5064318300	9,5 mm	R	300		SN5064111300	11,5 mm	R	300
SN5064318320	9,5 mm	R	320		SN5064111320	11,5 mm	R	320
SN5064318340	9,5 mm	R	340		SN5064111340	11,5 mm	R	340
SN5064318360	9,5 mm	R	360		SN5064111360	11,5 mm	R	360
SN5064318380	9,5 mm	R	380		SN5064111380	11,5 mm	R	380
SN5064318400	9,5 mm	R	400		SN5064111400	11,5 mm	R	400
SN5064318420	9,5 mm	R	420		SN5064111420	11,5 mm	R	420
SN5064318440	9,5 mm	R	440		SN5064111440	11,5 mm	R	440
SN5064328400	9,5 mm	L	300		SN5064221300	11,5 mm	L	300
SN5064328420	9,5 mm	L	320		SN5064221320	11,5 mm	L	320
SN5064328440	9,5 mm	L	340		SN5064221340	11,5 mm	L	340
SN5064328460	9,5 mm	L	360		SN5064221360	11,5 mm	L	360
SN5064328480	9,5 mm	L	380		SN5064221380	11,5 mm	L	380
SN5064328500	9,5 mm	L	400		SN5064221400	11,5 mm	L	400
SN5064328520	9,5 mm	L	420		SN5064221420	11,5 mm	L	420
SN5064328540	9,5 mm	L	440		SN5064221440	11,5 mm	L	440
SN5064118300	10 mm	R	300		SN5064113300	12 mm	R	300
SN5064118320	10 mm	R	320		SN5064113320	12 mm	R	320
SN5064118340	10 mm	R	340		SN5064113340	12 mm	R	340
SN5064118360	10 mm	R	360		SN5064113360	12 mm	R	360
SN5064118380	10 mm	R	380		SN5064113380	12 mm	R	380
SN5064118400	10 mm	R	400		SN5064113400	12 mm	R	400
SN5064110420	10 mm	R	420		SN5064113420	12 mm	R	420
SN5064110440	10 mm	R	440		SN5064113440	12 mm	R	440
SN5064220300	10 mm	L	300		SN5064223300	12 mm	L	300
SN5064220320	10 mm	L	320		SN5064223320	12 mm	L	320
SN5064220340	10 mm	L	340		SN5064223340	12 mm	L	340
SN5064220360	10 mm	L	360		SN5064223360	12 mm	L	360
SN5064220380	10 mm	L	380		SN5064223380	12 mm	L	380
SN5064220400	10 mm	L	400		SN5064223400	12 mm	L	400
SN5064220420	10 mm	L	420		SN5064223420	12 mm	L	420
SN5064220440	10 mm	L	440		SN5064223440	12 mm	L	440

Specifications

- Lateral Femoral Intramedullary Nail Instrument Set

Indications

- Lateral Femoral Nail with standard locking is indicated for fractures in the femoral shaft except subtrochanteric fractures femoral shaft in case of combination with femoral neck fractures: combined with double ipsilateral fractures indicated for fractures in the subtrochanteric section.

INTRAMEDULLARY CANNULATED LATERAL FEMUR NAIL

Lateral Femur Nail Cannulated Lag Screw

Code	Size	Length
SN5066230065	6,5	65
SN5066230070	6,5	70
SN5066230075	6,5	75
SN5066230080	6,5	80
SN5066230085	6,5	85
SN5066230090	6,5	90
SN5066230095	6,5	95
SN5066230100	6,5	100
SN5066230105	6,5	105
SN5066230110	6,5	110
SN5066230115	6,5	115
SN5066230120	6,5	120



Lateral Femur Nail Locking Screw

Code	Size	Length
SN5064205020	5 mm	20
SN5064205025	5 mm	25
SN5064205028	5 mm	28
SN5064205030	5 mm	30
SN5064205032	5 mm	32
SN5064205034	5 mm	34
SN5064205035	5 mm	35
SN5064205036	5 mm	36
SN5064205038	5 mm	38
SN5064205040	5 mm	40
SN5064205042	5 mm	42
SN5064205044	5 mm	44
SN5064205045	5 mm	45
SN5064205048	5 mm	48
SN5064205050	5 mm	50
SN5064205052	5 mm	52
SN5064205054	5 mm	54
SN5064205055	5 mm	55
SN5064205056	5 mm	56
SN5064205058	5 mm	58
SN5064205060	5 mm	60
SN5064205064	5 mm	64
SN5064205065	5 mm	65
SN5064205068	5 mm	68
SN5064205070	5 mm	70
SN5064205072	5 mm	72
SN5064205074	5 mm	74
SN5064205075	5 mm	75
SN5064205076	5 mm	76
SN5064205078	5 mm	78
SN5064205080	5 mm	80
SN5064205085	5 mm	85
SN5064205090	5 mm	90
SN5064205095	5 mm	95
SN5064205100	5 mm	100
SN5064205105	5 mm	105
SN5064205110	5 mm	110
SN5064205115	5 mm	115
SN5064205120	5 mm	120



Lateral Femur Nail Non-Cannulated Lag Screw

Code	Size	Length
SN5066240065	6,5	65
SN5066240070	6,5	70
SN5066240075	6,5	75
SN5066240080	6,5	80
SN5066240085	6,5	85
SN5066240090	6,5	90
SN5066240095	6,5	95
SN5066240100	6,5	100
SN5066240105	6,5	105
SN5066240110	6,5	110
SN5066240115	6,5	115
SN5066240120	6,5	120



Lateral Femur Nail End Cap

Code	Size
SN5007031000	0
SN5007031005	5 mm
SN5007031010	10 mm
SN5007031015	15 mm
SN5007031020	20 mm



TIBIAL NAIL CANNULATED

Code	Size	Length
SN5074208255	8,4 mm	255
SN5074208270	8,4 mm	270
SN5074208280	8,4 mm	285
SN5074208300	8,4 mm	300
SN5074208315	8,4 mm	315
SN5074208330	8,4 mm	330
SN5074208345	8,4 mm	345
SN5074208360	8,4 mm	360
SN5074208375	8,4 mm	375
SN5074208400	8,4 mm	400
SN5074209255	9 mm	255
SN5074209270	9 mm	270
SN5074209280	9 mm	285
SN5074209300	9 mm	300
SN5074209320	9 mm	315
SN5074209330	9 mm	330
SN5074209340	9 mm	345
SN5074209360	9 mm	360
SN5074209380	9 mm	375
SN5074209400	9 mm	400
SN5074210255	10 mm	255
SN5074210270	10 mm	270
SN5074210280	10 mm	285
SN5074210300	10 mm	300
SN5074210320	10 mm	315
SN5074210330	10 mm	330
SN5074210340	10 mm	345
SN5074210360	10 mm	360
SN5074210380	10 mm	375
SN5074210400	10 mm	400



Code	Size	Length
SN5074211255	11 mm	255
SN5074211270	11 mm	270
SN5074211280	11 mm	285
SN5074211300	11 mm	300
SN5074211320	11 mm	315
SN5074211330	11 mm	330
SN5074211340	11 mm	345
SN5074211360	11 mm	360
SN5074211380	11 mm	375
SN5074211400	11 mm	400
SN5074212270	12 mm	270
SN5074212280	12 mm	285
SN5074212300	12 mm	300
SN5074212320	12 mm	315
SN5074212330	12 mm	330
SN5074212340	12 mm	345
SN5074212360	12 mm	360
SN5074212380	12 mm	375
SN5074212400	12 mm	400

Specifications

- Tibial Intramedullary Nail Instrument Set

Indications

- Fractures in the tibial shaft as well as for metaphyseal and certain intraarticular fractures of the tibial head and the pilon tibial.

TIBIAL NAIL CANNULATED

Tibial Nail Locking Screw

Code	Size	Length
SN5076235025	4,8 mm	25
SN5076235030	4,8 mm	30
SN5076235035	4,8 mm	35
SN5076235040	4,8 mm	40
SN5076235045	4,8 mm	45
SN5076235050	4,8 mm	50
SN5076235055	4,8 mm	55
SN5076235060	4,8 mm	60
SN5076235065	4,8 mm	65
SN5076235070	4,8 mm	70
SN5076235075	4,8 mm	75
SN5076235080	4,8 mm	80
SN5076235085	4,8 mm	85
SN5076235090	4,8 mm	90
SN5076235095	4,8 mm	95
SN5076235010	4,8 mm	100



Tibial Nail Locking Screw

Code	Size	Length
SN5076234520	4,5 mm	20
SN5076234524	4,5 mm	24
SN5076234525	4,5 mm	25
SN5076234526	4,5 mm	26
SN5076234528	4,5 mm	28
SN5076234530	4,5 mm	30
SN5076234532	4,5 mm	32
SN5076234534	4,5 mm	34
SN5076234535	4,5 mm	35
SN5076234536	4,5 mm	36
SN5076234538	4,5 mm	38
SN5076234540	4,5 mm	40
SN5076234542	4,5 mm	42
SN5076234544	4,5 mm	44
SN5076234545	4,5 mm	45
SN5076234546	4,5 mm	46
SN5076234548	4,5 mm	48
SN5076234550	4,5 mm	50
SN5076234554	4,5 mm	54
SN5076234555	4,5 mm	55
SN5076234560	4,5 mm	60
SN5076234565	4,5 mm	65
SN5076234570	4,5 mm	70
SN5076234575	4,5 mm	75
SN5076234580	4,5 mm	80
SN5076234585	4,5 mm	85
SN5076234590	4,5 mm	90
SN5076234595	4,5 mm	95
SN5076234510	4,5 mm	100



Tibial Nail End Cap

Code	Size
SN5007032000	0
SN5007032005	5 mm
SN5007032010	10 mm
SN5007032015	15 mm
SN5007032020	20 mm



TIBIA ARTHRODESIS NAIL

Code	Size	Length
SN8501000150	10 mm	150
SN8501000175	10 mm	175
SN8501000200	10 mm	200
SN8501000225	10 mm	225
SN8501100150	11 mm	150
SN8501100175	11 mm	175
SN8501100200	11 mm	200
SN8501100225	11 mm	225



Arthrodesis Locking Screw

Code	Size	Length
SN5026230020	5.0 mm	20
SN5026230025	5.0 mm	25
SN5026230030	5.0 mm	30
SN5026230035	5.0 mm	35
SN5026230040	5.0 mm	40
SN5026230045	5.0 mm	45
SN5026230050	5.0 mm	50
SN5026230055	5.0 mm	55
SN5026230060	5.0 mm	60
SN5026230065	5.0 mm	65
SN5026230070	5.0 mm	70
SN5026230075	5.0 mm	75
SN5026230080	5.0 mm	80



Arthrodesis Nail End Cap

Code	Size
SN5026230085	0



Specifications

- Arthrodesis Tibial Nail Instrument Set

Indications

- Pertrochanteric fractures, Intertrochanteric fractures, High subtrochanteric fractures, Low subtrochanteric fractures, Femoral shaft fractures, Isolated or combined medial femoral neck fractures.

FLEXIBLE INTRAMEDULLARY NAIL

Code	Size	Length
SN5032202040	2,0 mm	400
SN5032202540	2,5 mm	400
SN5032203040	3,0 mm	400
SN5032203540	3,5 mm	400
SN5032204040	4,0 mm	400
SN5032204540	4,5 mm	400
SN5032205040	5,0 mm	400

Flexible Nail End Cap

Code	Size	Used with
SN5032202041	0	2&2.5 Ø Nails
SN5037030000	0	3-5 Ø Nails



Specifications

- Flexible Titanium Elastic Nail Instrument Set

Indications

- The management of diaphyseal and certain metaphyseal/-epiphyseal fractures of long bones in children and young adults.

09

ON THE HORIZON..

NEW PRODUCTS

DISTAL TIBIA ANTEROLATERAL LOCKING PLATE



PROXIMAL MEDIAL TIBIA LOCKING PLATE



ARTRODEZ LOCKING PLATE



3.5 VOLAR DISTAL RADIUS LOCKING PLATE



DISTAL POSTERIOR TIBIA LOCKING PLATE

DISTAL TIBIA PILON LOCKING PLATE

HERBER SCREWS SYSTEM



SANOMED

ADVANCED ORTHOPEDIC SOLUTIONS

SCREWS PLATES & NAILING

PRODUCT CATALOG 2025

SOLUTIONS

CONTACT SANOMED

- **WEB** www.sanomed.health
- **LOCATION** Koura, Lebanon
- **SCOPE** Orthopedic implants and surgical solutions

For professional use only. Product codes and technical details are presented as supplied by the manufacturer. Availability and regulatory status may vary. Contact SANOMED for current product and ordering information.

PRECISION. STRENGTH. HEALING.